

Text
MATERIALS LIST

| Mat. ID | CE | Material Name | Brand | Ingredients |
|---------|------|------------------------------------------|-----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 6 | 0123 | IPS E-MAX CERAM | IVOCLAR | SiO ₂ ,Al ₂ O ₃ ,ZnO ₂ ,Na ₂ O,K ₂ O,ZrO, CaO,P ₂ O ₅ , fluoride, pigments |
| 8 | 0123 | IPS E-MAX PRESS INGOTS | IVOCLAR | SiO ₂ ,Li ₂ O,K ₂ O,MgO,ZnO,Al ₂ O ₃ ,P ₂ O ₅ |
| 10 | 0120 | D-sign 98 | IVOCLAR | Au85.9%,Pt12.1%,Zn1.5%,In<1.0%,Ir<1.0% |
| 12 | 0120 | Lodestar | IVOCLAR | Au51.5%, Pd38.5%,In8.5%, Ga1.5%, Ru <1.0%, Re <1.0% |
| 13 | 0120 | Evolution Lite | IVOCLAR | Au40.3%,Pd39.3%,Ag9.2%,In9.3%,Ga1.8%,Ru<1.0% |
| 14 | 0120 | MiniGold | IVOCLAR | Au40.0%, Pd4.0%,Ag47.0%, Cu7.5%,Zn1.0%, In<1.0%,Ir<1.0% |
| 17 | 0197 | BEGO WIRONIT extrahart | BEGO | Co 63%,Cr 30%,Mo 5%,Si 1.1%,Mn,C |
| 18 | 1275 | Girobond® nbs | Girrbach | Co62.4%,Cr 25.5%,Mo 5.1%, W5.2%,Si1.1%,Nb, Fe,N<1% |
| 19 | 0120 | Self Curing-Pink | Vertex | Polymethylmethacrylate, plasticiser, pigments and catalysts, 0.5-1.5%Benzoylperoxide |
| 20 | 0483 | Orthocryl Powder-Clear | Dentaurum | Methylmethacrylate95%,Butandiole,Dimethacrylate<5% |
| 22 | 0123 | Ivostar | IVOCLAR | Methylmethacrylate |
| 23 | 0483 | Ropak uv-fluid | BREDDENT | |
| 25 | 0044 | M.L. Primer | SHOFU | Azeton80%,other |
| 26 | 0470 | Preci-Vertex | CEKA | |
| 27 | 0483 | FINO AUTOBAOND | FINO | ATT Adhesive |
| 28 | 0470 | VALPLAST | VALPLAST | Nylon |
| 29 | 0120 | HIGH FUSING WHITE CERAMIC SOLDER | IVOCLAR | Au45.0%,Pd12.4%,Ag41.5%, In1.0%, Li<1.0%,Ru<1.0% |
| 30 | 0120 | SPECIAL HIGH FUSING WHITE CERAMIC SOLDER | IVOCLAR | Au47.0%, Pd10.3%,Ag41.0%, B<1.0%, In1.4%, Ru<1.0%, Other<1.0% |
| 31 | 0120 | UNIVERSAL SOLER 1015W | IVOCLAR | Au18.5%, Pd5.99%, Ag72.5%, In3.0%, Ir<1.0% |
| 32 | 0120 | UNIVERSAL SOLDER PKF | IVOCLAR | Au48.8%, Ag40.5%, Zn7.3%, Pt2.8%, Other<1.0% |
| 38 | 0123 | E-MAX GLAZE | IVOCLAR | Ethandiol |
| 39 | 0123 | E-MAX CERAM MARGIN | IVOCLAR | SiO ₂ ,Al ₂ O ₃ ,ZnO ₂ ,Na ₂ O,K ₂ O,ZrO, CaO,P ₂ O ₅ , fluoride, pigments |

| Mat. ID | CE | Material Name | Brand | Ingredients |
|---------|-----------|---------------------------|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 53 | 0483 | MK1 | DENTAL-ATTACHMEN | |
| 55 | 0044 | TK-SNAP | SI-TEC | |
| 56 | 0123 | E-MAX CERAM GINGIVA | IVOCLAR | SiO ₂ ,Al ₂ O ₃ ,ZnO ₂ ,Na ₂ O,K ₂ O,ZrO, CaO,P ₂ O ₅ ,fluoride,pigments |
| 58 | 0124 | TITANIUM CAD CAM ABUTMENT | GEO | Titanium80-100%,Molybdenum1-15%, Aluminum<10%, Zirconium<10%,Iron<5%,Chromium<5%,Copper<5%,Tin<5%,Cobalt0-Trace,Manganese0-Trace, Nickel0-Trace,Tantalum0-Trace, Tungsten0-Trace,Vanadium0-10% |
| 60 | 1251 | BALL ATTACHMENT | CEKA | |
| 63 | 0124 | TITANIUM CERAMIC | VITA | SiO ₂ 63-68%,Na ₂ O 5-7%,Cao 1.5-2.5%,K ₂ O 8-10%,Al ₂ O ₃ 8-10%,BaO 0.3-1%,ZrO ₂ <0.5%,SnO ₂ <0.5%,B ₂ O ₃ 7-9% |
| 64 | 0120 | PMMA DISK | YAMAHACHI | |
| 66 | 0483 | DENTAURUM CLASP | DENTAURUM | NICKEL FREE |
| 71 | 0483 | BALL RETAINER CLASPS | DENTAURUM | |
| 72 | 0197 | Upzir solid zirconia | UPCERA | ZrO ₂ +HfO ₂ +Y ₂ O ₃ ≥ 99%,Y ₂ O ₃ =4.5~6.0%,HfO ₂ ≤ 5%,Al ₂ O ₃ ≤ 0.5%,other ≤ 0.5% |
| 74 | | Resilient Acrylic | Bosworth Trusoft | |
| 77 | 0297 | mega TRIM Provi Implant | megadental | |
| 78 | 0483 | Remanium Draht/Wire | Dentaurum | |
| 79 | 0483 | Expansion Screw | Dentaurum | |
| 80 | 0344 | CZR | NORITAKE | |
| 81 | 0120 | Harmony Hard | IVOCLAR | Au74.0%,Pd3.8%,Ag12.0%,Cu9.0%,Zn<1.0%, In<1.0%, Ir<1.0% |
| 82 | 0120 | SPARTAN PLUS | IVOCLAR | Au2.0, Pd78.8, Cu10.0, Ga9.0, Ir<1.0, Ge<1.0, Li<1.0 |
| 83 | 93/42/EWG | Silensor-sl Teilekarten | ERKODENT | polyethylene |
| 85 | 93/42/EWG | ERKODUR | ERKODENT | Ethylene terephthalate,Ethylene glycol,Vinyl acetate, Ethylene polymerization |
| 86 | 93/42/EWG | ERKOFLEX-BLEACH | ERKODENT | Ethylene terephthalate,Ethylene glycol,Vinyl acetate, Ethylene polymerization |
| 87 | 0482 | DDcubeX2 | Dental Direkt materials | ZrO ₂ +HfO ₂ +Y ₂ O ₃ ≥ 99%; Y ₂ O ₃ < 10%; Al ₂ O ₃ ~ 0.05%; other ≤ 0.05% |
| 88 | 289250 | Vario Ball Snap | BREDDENT | |
| 90 | 0123 | IPS Classic | IVOCLAR | SiO ₂ 40~65wt.%, Al ₂ O ₃ ,B ₂ O ₃ , BaO, CaO, CeO ₂ , K ₂ O,MgO, Na ₂ O, P ₂ O ₅ , SnO ₂ , TiO ₂ , ZrO ₂ and pigments |
| 91 | 0123 | SR Nexco Opaquer | IVOCLAR | IVOCLAR VIVADENT |

| Mat. ID | CE | Material Name | Brand | Ingredients |
|---------|-----------|-------------------------------|---------|----------------------------------------|
| 92 | 0123 | SR Nexco Paste | IVOCLAR | IVOCLAR VIVADENT |
| 95 | 0123 | E-MAX GLAZE | IVOCLAR | Ethandiol |
| 96 | 0123 | IPS E-MAX CERAM GLAZE SPRAY | IVOCLAR | Ethandiol |
| 97 | 0124 | Reinforcement mesh nets | RENFERT | |
| 98 | 0120 | THERMOSENS RIGID | VERTEX | Thermoplastic Dental Resin |
| 102 | 524151MR5 | GIALLOY TITANIUM ALLOY BLANKS | SRL | Ti 90.04%, Al5.94%, V4.0%,Fe0.02%,O 0% |